OVERVIEW

TECHNICAL SPECIFICATION

MANUAL/LINKS

The Adjustable Bath & Shower Seat can be used in a shower stall or bathtub. The seat is constructed of durable plastic with easy-to-grip handles and a notch to hold a hand-held shower spray. The aluminum legs are flared to help prevent tipping, and the slip-resistant rubber tips provide added stability.

- Comfortable contoured seat and backrest. The bath chair combines shower safety and comfort with a contoured design for long periods of sitting.
- Holds up to 350 lbs. The heavy-duty shower chair is made of durable molded plastic with aluminum legs. The aluminum frame flares help distribute weight for safety and stability.
- Built-in handles provide stability. The shower seat for seniors features side handles for a grippable support surface. This makes it an ideal handicap shower seat by preventing falls and slips.
- Height adjustable from 14" to 20". The adjustable shower chair easily adjusted with the push-of-a-button on each leg. Having the correct shower seat height is essential.
- Easy, tool-free assembly. The shower chair for elderly assembles in minutes with no tools needed. It can be taken apart to be used as a convenient travel shower chair.
- Available with or without back. Find the bath and shower chair with a back here.

TECHNICAL SPECIFICATION

MANUAL/LINKS

Weight Capacity: 350 lbs.

Height Adjustable: Yes

Back: No

Handles: No

Package Dims (h x w x d): 22" x 13.5" x 4"

Product Dims (h x w x d): 13.5" - 20.5" x 20" x 18"

Product Weight: 4.1 lbs.

Package Weight: 6.3 lbs.